

# Microbial Preservation Plate **96Well**

Under normal temperature and pressure  
**Our record preservation time  
is 5 years**

\*Preservation time varies depending on sample type.



- **Dramatic savings per chip!**
- **Made from Polypropylene - it is easy to cut !**
- **Convenient well layout to match  
your 96 well plate!**

## Features:

Microbial strains such as *E.coli* and yeast can be dry-preserved on a paper chip coated with a preservation agent.

Samples can be sealed individually with an accessory protection seal to prevent contamination and inactivation.

This enables compact and efficient preservation at normal room temperature - it is no longer necessary to use deep freeze or dry ice for transportation.

Note : When handling the specific pathogens, please act in accordance with the applicable laws and regulations in your country.



**Preservable in file folders,  
saving you storage space!**

Cat. No.	Item	Unit
176-531S	Microbial preservation plate for Bacteria, 96well	<small>Sterilized</small> 5 plates / bag
176-551S	Microbial preservation plate for Fungi, 96well	<small>Sterilized</small> 5 plates / bag

\* Specification is subject to change without notice.

# 3 Features



## 1 Cost reduction!

Dramatic increase in the number of samples that can be preserved per plate.



## 2 Easy to cut as you like!

The plate is made from soft PP material that can be easily cut to any required size.


It is very convenient when you have more than one type of sample and you need to prevent contamination.

Watson preservation plates have a wide range of applications.




Easy to cut!

### Preservation Method




- ① Let the paper chip absorb a liquid sample.  
A: For a liquid sample  
B: For a sample from agar
- To the level  $5\mu\text{L}$  of solution volume per well

- ② Dry the sample. (Dry for at least 60 mins. Reduced pressure drying is recommended.)  
③ Seal the plate with the protection seal and store it in room temperature.



### Recovery Method



- ① Peel off the protection seal.  
② Place the paper chip into a container.  
③ In the case of agar, rub the paper chip onto agar.

# Q & A

**Q.** Correlation between preservation time and storage temperature?

**A.** Please store at the lower temperature for longer preservation. Some samples can be stored at the normal temperature if it is for short period. Please see our lab data at the back page for your reference.

**Q.** What is the shelf life of the products before unsealing and use?

**A.** 1 year after the manufacturing date.

**Q.** Can I store the unsealed product?

**A.** Please do not store the unsealed product without using it. Sterile condition is lost once it is unsealed. Please unseal the products in a clean bench or a safety cabinet just before the use.

**Q.** Is it better to store the plates at -80°C rather than -20°C ?

**A.** Please avoid storing the products at lower than -20°C. The protection seal may not function under the temperature lower than -20°C.

**Q.** What is the recommended method for storing plates ?

**A.** Bundle 2 plates with a rubber band so that protection seals face to each other and place it in a ziplock bag. This prevents the protection seal from rubbing and getting damaged by another plate.

## ■ Microbial Preservation Plate Product Line



1 Well



3 Well

Cat. No.	Item	Unit
176-231S	Microbial preservation plate for Bacteria, 1well	Sterilized 5 plates / unit
176-251S	Microbial preservation plate for Fungi, 1well	Sterilized 5 plates / unit

Cat. No.	Item	Unit
176-331S	Microbial preservation plate for Bacteria, 3well	Sterilized 5 plates / unit
176-351S	Microbial preservation plate for Fungi, 3well	Sterilized 5 plates / unit

For preservation of nucleic acid (DNA, RNA) or blood samples....

### Preservation Plate (PVP)

\* Nylon paper chip type

\*Cellulose paper chip type



Preservation Plate has been developed from the study result of MEXT's Intellectual Cluster Formation Project <Tokushima Region Noji group (The University of Tokushima)>. This project is supported by METI's New Cooperation Measure Subsidy.

- Our in-house verification experiments have been implemented in aerobic condition only.  
For preservation of anaerobic bacteria, please conduct pretest in advance.

#### ● Recommended dry time

(For ambient condition of 23°C, Humidity 50%)

[96 well] 60 minutes or more  
[1/3 well] 90 minutes or more



### Precautions

#### ● When handling certain pathogenic samples, please handle it in compliance with pertinent laws.

- Do not use Preservation Plate for other purpose than study.
- Wear gloves and a mask when you operate.
- Preservation Plate is disposable. Do not reuse.
- Do not autoclave.
- Store Preservation Plate avoiding high temperature and humidity after unsealing.

- Please avoid light, dusts, high humidity for sample storage space.
- Preservation time may vary depending on purity and/or storage condition of a sample.
- Conduct half-life test to grasp preservation time.  
[Half-life : t(1/2) = ln2 / {ln(100) - ln(Survival rate after a month)}]

