

WATSON®

Dispo Pipette

disposable pipette (5ml/10ml)

Experience high precision graduation achieved by our injection molding technology



Readability of the graduation enhances as aspirating a solution as the contrast increases.



Cat.No.	Item	Packaging
1650-005S	5ml dispo pipette	Sterilized 150 pcs.(individual packaging)
1650-010S	10ml dispo pipette	Sterilized 100 pcs.(individual packaging)

watsonbiolab.com

Short scale,
for the ease of
use inside
a clean bench

Body: PS / Adapter:PP
DNase/RNase free
Pyrogen test results:
Endotoxin level less
than 0.05EU/ml
*according to our internal test

For Research

Made in JAPAN

Dispo Pipette



Something Different.
WATSON BIO LAB
MADE IN JAPAN SINCE 1988

Manufactured and sold by:

FUKAEKASEI CO., LTD.

Head Office: 2-2-7 Murotani, Nishi-ku, Kobe


651-2241, JAPAN

TEL +81-78-991-4477

FAX +81-78-991-4491

E-mail: info@watson.co.jp

<https://watsonbiolab.com>



※The specification and the price etc. of the product might change without a previous notice.